

> **Welcome to the Revolution**

The Exprt™ Precision System–Revision Knee redefines total knee arthroplasty, providing advanced technology that translates into greater efficiency in the operating suite. The system’s streamlined, compact design reduces turn times, minimizes waste, and has proven implant design technology - all at 40-70% the cost of comparable knee revision systems.

> **EXPRT Features:**

Built in Augments

4mm distal and posterior augments built into femoral component compensate for bone loss

Modular Stem Extensions (Straight)

75mm straight and 65mm tapered stems may be used interchangeably on the femur and tibia

Single Radius Design

4 asymmetric femoral sizes with anatomic single radius to provide stability and flexibility

Tibial Insert

12mm, 14.5mm, 17mm and 19.5mm tibial insert thicknesses

Reduced wear

e+ poly is formulated to reduce long-term wear and neutralizes the negative effects of free radicals

Full Tibial Augment

10mm optional augment compensates for bone loss

Thicker Baseplates

4 tibial sizes and 4 symmetric 6mm tibial baseplates



Modular Stem Extensions (Tapered)

Threaded stem extensions for a stable construct



> Value Driven Innovation

The Exprt™ Revision Knee features two instrument trays vs. competitors' eight trays, saving valuable time and space in the operating room.



> Instrumentation:

- Only 2 Instrument Trays for the back-table
- Same Distal/ Tibial Cutting Guide
- One self-centering reamer for both femoral and tibial preparation

